



<p><b>Transfer Board, Curved</b></p> <p>Plastic curved 'Transfer Board', light weight, easy to manage with two anti-slip pads for added security.</p> <ul style="list-style-type: none"><li>• <b>Size : 75 cm x 25 cm (29" x 10")</b></li><li>• <b>Must have safe working load of: 285 kg (44 stone)</b></li></ul> <p><b>Price band A</b></p> <p style="text-align: right;">Aug 11</p>	<p><b>Order Code: MOV 0503</b></p> 
<p><b>Transfer Board, Straight</b></p> <p>Slick surface reduces energy expense of user and / or effort of carer. Lightweight construction of durable polyethylene for flexibility and maximum surface contact, perfect for transfers with different elevations. Friction tape on the bottom keeps it safely in place. Hand cut outs allow for easier positioning and removal once transfer is complete.</p> <ul style="list-style-type: none"><li>• <b>Size : 76cm x 32cm (30" x 12.5")</b></li><li>• <b>Must have safe working load of : 190 kg (29 stone)</b></li></ul> <p><b>Price band B</b></p> <p style="text-align: right;">Aug 11</p>	<p><b>Order Code: MOV 0504</b></p>  <p>Also supplied in yellow and green</p>

### Swivel Seat / Transfer Disc

Features include: two layer low friction disc; slip-resistant base; Soft, washable top. Suitable for pivot transfer.

Gives greater mobility when transferring to and from cars. Comes complete with car seat belt loop

- Overall Diameter approx 340mm
- Height max 50mm



**When using in cars always ensure retaining strap is used and securely fitted to car seat back.**

**Price band A**

Aug 11



**This item must be ordered using a retail prescription**



**Retail Model Order Code: MO01**

**Retail Model** - see Devon Tariff document for details.

**CES Order Code: MOV 0614**

### Transfer Belt

Soft padded handling belt with a non-slip lining to help prevent the belt sliding upwards during transfer (sit to stand or assisted walking)

**Small : 56cm to 81cm (22" to 31"), order code : MOV 0505**

**Medium : 76cm to 117cm (30" to 46"), order code : MOV 0506**

**Large : 120cm to 142cm (47¼" to 56"), order code : MOV 0507**



**Not to be used for lifting patients**

**Price band A**

Jun 08

### Order Codes :

**See detail**



### Turntable

Patient turner to move people into a standing position from a seated position or vice versa, or to transfer from one seated position to another. The turner has adjustable, padded lower leg supports and adjustable handle height.

- **Must have safe working load of: 150 kg (24 stone)**

Price band D

Aug 11

CES Order Code MOV 0627



### Return 7500 Transfer Device

The Return 7500 enables a pivoting transfer combined with the facility to transport service users short distances.



**Should be issued by Physiotherapist, O.T. or Nurse, trained in moving and handling, following an assessment.**

A Return handling belt is recommended for use with this product but must be ordered separately.

- **Low instep fits between wheelchair wheels.**
- **Adjustable knee brace pads.**
- **Must have safe working load of: 150 kgs (24st)**

Price band D

Aug 11

CES Order Code: MOV 0643



[See Prescriber guidelines](#)

### Return Handling Belt

A handling belt equipped with several strong handles to guide the user.

User can be safely strapped and secured to the stem of the Return 7500 transfer device.

Available in various sizes: Size corresponds to user's hip measurement.

- Small: 850-1050mm
- Med: 1050 – 1250mm
- Large: 1200 – 1450mm

**Price band A**

Jan 11



**CES Order Codes :**

**Small – MOV0749**

**Medium – MOV0655**

**Large – MOV0750.**

### Transfer Sheet, Multi Glide

The flat slide sheet is an alternative to the tubular slide sheets. It has handles on 3 sides of the slide and has the same friction resistance as the tubular slide sheets. It facilitates transfer and positioning of the user in multiple directions.

- **Size 200cm (79") x 71cm (28")**

**PRICE BAND A**

Jun 08

**CES Order Code : MOV 0512**



**Transfer Sheet, Multi Glide**

Low friction nylon fabric roller. It facilitates transfer and positioning of the user in multiple directions.

For sizes and order codes:

[see chart](#)








**Price band B**

Order codes: [see chart](#)

Jun 08



### Slide Sheet Sizing Chart

Direction of slide sheet	Measurements	Code
	72cm x 71cm	<b>MOV 0510</b>
	122cm x 71cm	<b>MOV 0559</b>
	122cm x 100cm	<b>MOV0509</b>
	145cm x 71cm	<b>MOV0569</b>
	200cm x 70cm	<b>MOV 0564</b>
	200cm x 140cm	<b>MOV 0508</b>
	200cm x 71cm	<b>MOV 0511</b>

## Wendylett System

A system of in situ slide and draw sheets which work in combination to facilitate positioning and transfers.



**Should be issued by Physiotherapist, O.T. or Nurse, trained in moving and handling, following an assessment.**

The following items are available

- **Single Base sheet**
- **Fitted full length (maxi)**
- **Flat three quarter length (midi)**
- **2 Way draw sheet**
- **4 Way draw sheet**

[See prescriber guidelines](#)

Price band D

Jan 11



**CES Order Codes:**

**Fitted single base sheet : MOV0757**

**Flat midi base sheet: MOV0886**

**2 Way draw sheet: MOV0887**

**4 Way draw sheet: MOV0608**