

Prescribing Guidelines for Return7500

Where MOV 0627 doesn't work

Indications for use	Indications for non use	Further action to be taken	Alternatives
<ul style="list-style-type: none"> • To be issued by nurse, O.T or physiotherapist trained in moving and handling following an assessment • User able to stand but unable to step around • User needs help to pull himself to standing. • Difficulty in getting furniture in the correct position e.g chair. • Possibility of reducing number of care assistants required. • Handling belt attached to Return enables user to safely be moved in standing. 	<ul style="list-style-type: none"> • User can stand and step around. • User non weight bearing • MOV 0627 type turner would be appropriate, i.e. user can stand while pulling on grab handle and can support themselves by holding the grab rail during manoeuvre through 90 to other furniture. • Movement of Return destabilises user too much • Some floor surfaces e.g. thick carpets make Return difficult to move. • Return can not be moved near enough to seating surface for safe stand to sit. 	<ul style="list-style-type: none"> • Thorough moving and handling assessment by O.T. Physiotherapist or Nurse. • Assess whether 'handling belt' required to go with Return and correct size. • Check all care assistants can use equipment correctly/safely. 	<ul style="list-style-type: none"> • MOV 0627 turntable • Side transfer • Transfer board • Stand assist hoist • Macnish Chair raiser • Hoist